

Native Staphylococcus aureus Protein A (Alkaline Phosphatase) ab7345

描述

产品名称	Native金黄色葡萄球菌蛋白A (Alkaline Phosphatase)
纯度	> 95 % Protein A purified.
表达系统	Native
Accession	<u>P38507</u>
蛋白长度	Full length protein
无动物成分	No
性质	Native
种属	Staphylococcus aureus
预测分子量	56 kDa
偶联物	Alkaline Phosphatase

技术指标

Our **Abpromise guarantee** covers the use of **ab7345** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot Dot blot ELISA Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) Immunohistochemistry (Frozen sections)
形式	Liquid
补充说明	Protein A Concentration: 1.0 mg/ml (by UV absorbance at 280 nm) This product has been assayed against 1.0 ug of Human IgG in a standard capture ELISA using pNPP (p-nitrophenyl phosphate) as a substrate for 30 minutes at room temperature. A working dilution of 1:500 to 1:2,000 of the stock concentration is suggested for this product. Optimal titers for other applications should be determined by the researcher.

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C. Do Not Freeze.

Preservative: 0.1% Sodium azide

Constituents: 0.00136% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl, 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, 1.5% BSA

BSA is IgG and Protease free

常规信息

相关性

Protein A is a surface protein of *S.aureus* which binds IgG molecules by their Fc region. In serum, the bacteria will bind IgG molecules in the wrong orientation on their surface which hinders opsonization and phagocytosis. Mutants of *S. aureus* lacking protein A are more efficiently phagocytosed in vitro, and mutants in infection models have diminished virulence. Due to its affinity for the Fc region of many mammalian immunoglobulins, protein A is considered a universal reagent in biochemistry and immunology.

细胞定位

Secreted; cell wall.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors