

NLRC5 peptide ab196047

描述

产品名称	NLRC5多肽
无动物成分	No
性质	Synthetic

技术指标

Our **Abpromise guarantee** covers the use of **ab196047** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Blocking - Blocking peptide for Anti-NLRC5 antibody (<u>ab105411</u>)
形式	Liquid
补充说明	This is the blocking peptide for <u>ab105411</u>

制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. pH: 7.20 Preservative: 0.02% Sodium azide Constituents: PBS, 0.1% BSA
--------	--

常规信息

功能	Probable regulator of the NF-kappa-B and type I interferon signaling pathways. May also regulate the type II interferon signaling pathway. Plays a role in homeostatic control of innate immunity and in antiviral defense mechanisms.
组织特异性	Expressed in spleen, thymus, lung, brain, tonsil, heart and prostate.
序列相似性	Belongs to the NLRP family. Contains 26 LRR (leucine-rich) repeats. Contains 1 NACHT domain.
细胞定位	Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors