

# Prostate Specific Antigen peptide ab41421

## 1 References

### 描述

产品名称	Prostate Specific Antigen多肽
纯度	> 95 % SDS-PAGE.
无动物成分	No
性质	Synthetic

### 技术指标

Our **Abpromise guarantee** covers the use of **ab41421** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	Western blot
	Blocking
形式	Liquid

### 制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
	pH: 6.75
	Constituents: 0.001% Tween, 0.714% HEPES, 0.0584% EDTA, 0.87% Sodium chloride

### 常规信息

功能	Hydrolyzes semenogelin-1 thus leading to the liquefaction of the seminal coagulum.
序列相似性	Belongs to the peptidase S1 family. Kallikrein subfamily. Contains 1 peptidase S1 domain.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors