

## Product datasheet

# Recombinant BK polyomavirus, strain AS, Major Capsid VP1 protein ab74567

[2 References](#) [1 图像](#)

### 描述

产品名称	重组BK polyomavirus, strain AS, Major Capsid VP1蛋白
纯度	> 90 % SDS-PAGE. Purified by ultracentrifugation.
表达系统	Saccharomyces cerevisiae
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	BK virus
序列	<p>MAPTKRKGECPGAAPKKPKPEVQVPKLLIKGGVEVLEVKTV            DAITEVEC            FLNPEMGDPDDNLRGYSQHLSAENAFESDSPDRKMLPCYSTA            RIPLPNLN            EDLTCGNLLMWEAVTVKTEVIGITSMLNLHAGSQKVHENG            GGGKPVQGSNF            HFFAVGGDPLEMQGVLNMYRTKYPQGTITPKNPTAQSVMNT            DHKAYLDK            NNAYPVECIWIPDPSRNTYFGTYTGGENVPPVLHVTNTAT            TVLLDEQG            VGPLCKADSLYVSAADICGLFTNSSGTQQRGLARYFKIRLR            KRSVKNPY            PISFLLSDLINRRTQKVDGQPMYGMESQVEEVRVFDGTEQLP            GDPDMIRY IDRQGQLQTKMV</p>

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab74567** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	Western blot
	SDS-PAGE

**形式** Lyophilized

## 制备和贮存

---

**稳定性和存储** Shipped at 4°C. Store at 4°C prior to reconstitution. Store at -20°C long term. Reconstitute for long term storage.

Constituent: PBS

**复溶** Reconstitute with deionized H<sub>2</sub>O. With the addition of glycerol store at -20°C. After reconstitution can only be stored at 4°C for 2-3 days.

## 常规信息


---

**相关性** BK polyomavirus was first isolated in 1971 in the urine of a renal transplant patient, initials B.K. BK polyomavirus, Major Capsid VP1 forms the icosahedric capsid of the virus protecting the viral genome. This capsid is composed of 360 VP1, assembled in 72 pentamers. Each pentamer are linked via the C-terminal end of VP1.

**细胞定位** Virion. Nucleus

## 图片

---



SDS-PAGE showing ab74567 at approximately 44kDa (3.2µg/lane).

Lane 1 shows the molecular weight ladder. From the bottom: 14.4, 18.4, 25.0, 35.0, 45.0, 66.2, 116.0 kDa.

SDS-PAGE - Recombinant BK polyomavirus, strain AS, Major Capsid VP1 protein (ab74567)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Promise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors