

Recombinant Cytomegalovirus Glycoprotein B ab43040

4 References

描述	
产品名称	重组Cytomegalovirus Glycoprotein B
生物活性	Immunoreactive with sera of CMV infected individuals.
纯度	> 95 % SDS-PAGE. Sephadex derived purification.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
额外的序列信息	The protein contains the Cytomegalovirus (CMV) gB immunodominant regions.
描述	重组Cytomegalovirus Glycoprotein B

技术指标	
Our Abpromise guarantee covers the use of ab43040 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	ELISA Western blot
形式	Liquid

制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.20 Constituents: 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol

常规信息	
相关性	Human cytomegalovirus (HCMV) glycoprotein B (gB) is an abundant virion envelope protein that has been shown to be essential for the infectivity of HCMV. HCMV gB is also one of the most immunogenic virus-encoded proteins, and a significant fraction of virus neutralizing antibodies are

directed at gB.

细胞定位

Virus envelope.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors