

# Recombinant Cytomegalovirus UL80a protein ab43042

## 描述

产品名称	重组Cytomegalovirus UL80a蛋白
生物活性	Immunoreactive with sera of CMV infected individuals.
纯度	> 95 % SDS-PAGE. Sephadex derived purification.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
额外的序列信息	The protein contains the CMV Pp38 (UL80a) immunodominant regions.

## 技术指标

Our **Abpromise guarantee** covers the use of **ab43042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA Western blot
形式	Liquid

## 制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. pH: 7.20 Constituents: 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol
--------	--

## 常规信息

相关性	Assembly of many spherical virus capsids is guided by an internal scaffolding protein or group of proteins that are often cleaved and eliminated in connection with maturation and incorporation of the genome. In cytomegalovirus one of the proteins that contributes to this scaffolding function is the maturational protease precursor (pUL80a).
-----	---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors