

Recombinant Dog GRP94 protein ab92290

[2 References](#) [2 图像](#)

描述

产品名称	重组狗GRP94蛋白
纯度	> 90 % SDS-PAGE. > 90% pure as determined by SDS-PAGE and Western Blot analysis
表达系统	Baculovirus infected Sf9 cells
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Dog

技术指标

Our **Abpromise guarantee** covers the use of **ab92290** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Western blot
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: 0.418% MOPS, 0.0146% EDTA
--------	--

常规信息

功能	Molecular chaperone that functions in the processing and transport of secreted proteins. Functions in endoplasmic reticulum associated degradation (ERAD). Has ATPase activity.
序列相似性	Belongs to the heat shock protein 90 family.
细胞定位	Endoplasmic reticulum lumen. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

图片



SDS-PAGE - Recombinant Dog GRP94 protein
(ab92290)

SDS-PAGE Analysis of GRP94 protein (ab92290) detected by
Coomassie Stain:


Lane 1: Mwt marker

Lane 2: 0.5 μ g

Lane 3: 1 μ g

Lane 4: 2 μ g

Lane 5: 5 μ g



Western blot - Recombinant Dog GRP94 protein
(ab92290)

Recombinant Dog GRP94 protein (ab92290) at 0.1 μ g

Predicted band size: 92 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors