

Recombinant Dog Prorenin protein ab92555

描述

| | |
|-------|--|
| 产品名称 | 重组狗Prorenin蛋白 |
| 纯度 | > 95 % SDS-PAGE. Purified by chelated metal affinity chromatography |
| 表达系统 | HEK 293 cells |
| 蛋白长度 | Full length protein |
| 无动物成分 | No |
| 性质 | Recombinant |
| 种属 | Dog |

技术指标

Our **Abpromise guarantee** covers the use of **ab92555** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----|----------|
| 应用 | SDS-PAGE |
| 形式 | Liquid |

制备和贮存

| | |
|--------|--|
| 稳定性和存储 | Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituent: 0.79% Tris HCl |
|--------|--|

常规信息

| | |
|------|---|
| 功能 | Renin is a highly specific endopeptidase, whose only known function is to generate angiotensin I from angiotensinogen in the plasma, initiating a cascade of reactions that produce an elevation of blood pressure and increased sodium retention by the kidney. |
| 疾病相关 | Defects in REN are a cause of renal tubular dysgenesis (RTD) [MIM:267430]. RTD is an autosomal recessive severe disorder of renal tubular development characterized by persistent fetal anuria and perinatal death, probably due to pulmonary hypoplasia from early-onset oligohydramnios (the Potter phenotype). Defects in REN are the cause of familial juvenile hyperuricemic nephropathy type 2 (HNFJ2) |

[MIM:613092]. It is a renal disease characterized by juvenile onset of hyperuricemia, slowly progressive renal failure and anemia.

序列相似性

Belongs to the peptidase A1 family.

细胞定位

Secreted. Membrane. Associated to membranes via binding to ATP6AP2.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors