

Recombinant E. coli Disulfide-bond isomerase protein ab51284


1 图像

描述	
产品名称	重组E. coli Disulfide-bond isomerase蛋白
纯度	> 95 % SDS-PAGE. This protein was purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Escherichia coli
序列	DDAAIQQTLA KMGIKSSDIQ PAPVAGMKTV LTNSGVLYIT DDGKHIIQGP MYDVSGTAPV NVTNKMLLKQ LNALEKEMIV YKAPQEKHVI TVFTDITCGY CHKLHEQMAD YNALGITVRY LAFPRQGLDS DAEKEMKAIW CAKDKNKAFD DVMAGKSVAP ASCDVDIADH YVLGVQLGVS GTPAVVLSNG TLVPGYQPPK EMKEFLDEHQ KMTSGK
氨基酸	21 to 236
描述	重组 <i>E. coli</i> Disulfide-bond isomerase蛋白
技术指标	
Our Abpromise guarantee covers the use of ab51284 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.242% Tris, 0.0584% EDTA

常规信息

相关性	Dsb proteins (DsbA, DsbB, DsbC, and DsbD) catalyze formation and isomerization of protein disulfide bonds in the periplasm of <i>Escherichia coli</i> . Disulphide bond isomerase (DsbC), a member of the thioredoxin superfamily, acts as a disulfide isomerase and converts aberrant disulfide bonds to correct ones. It also acts by transferring its disulfide bond to other proteins and is reduced in the process. DsbC is reoxidized by a yet uncharacterized protein.
细胞定位	Periplasm

图片



ab51284 run on a 15% SDS-PAGE gel with molecular weight markers.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors