

Recombinant E. coli Ecotin protein ab103303

1 图像

描述

产品名称	重组E. coli Ecotin蛋白
纯度	> 95 % SDS-PAGE. ab103303 was purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>P23827</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Escherichia coli
序列	MGSSHHHHHHSSGLVPRGSHMAESVQPLEKIAPYPQAEKG MKRQVIQLTP QEDESTLKVELLIGQTLEVDCNLHRLGGKLENKTLEGWGYDY YVFDKVSS PVSTMMACPDGKKEKKFVTAYLGDAGMLRYNSKLPVVYTPD NVDVKYRV WKAEEKIDNAVVR
预测分子量	18 kDa including tags
氨基酸	21 to 162
标签	His tag N-Terminus
描述	重组 <i>E. coli</i> Ecotin蛋白

技术指标

Our **Abpromise guarantee** covers the use of **ab103303** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
补充说明	Previously labelled as E. coli serine protease inhibitor.

制备和贮存

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.29% Sodium chloride

常规信息


相关性

E. coli serine protease inhibitor (Ecotin) inhibits pancreatic serine proteases: inhibits chymotrypsin, trypsin, elastases, factor X, kallikrein as well as a variety of other proteases. The strength of inhibition does not appear to be correlated with a particular protease specificity.

细胞定位

Periplasm.

图片



15% SDS-PAGE analysis of 3µg ab103303.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors