

### Recombinant *E. coli* fimC protein ab104677

#### 1 图像

#### 描述

<b>产品名称</b>	重组 <i>E. coli</i> fimC蛋白
<b>纯度</b>	> 95 % SDS-PAGE. ab104677 is purified using conventional chromatography techniques.
<b>表达系统</b>	<i>Escherichia coli</i>
<b>Accession</b>	<b><u>P31697</u></b>
<b>蛋白长度</b>	Full length protein
<b>无动物成分</b>	No
<b>性质</b>	Recombinant
<b>种属</b>	<i>Escherichia coli</i>
<b>序列</b>	<pre> MGSSHHHHHSSGLVPRGSHMGVALGATRVIYPAGQKQEQLA VTNNDENS TYLIQSWVENADGVKDRFIVTPPLFAMKGGKENTLRILDAT NNQLPQDR ESLFWMNVKAIPSMDSKLTENTLQLAIISRIKLYYRPAKLA LPPDQAAE KLRFRRSANSLTLINPTPYLLTVTELNAGTRVLENALVPPMG ESTVKLPS DAGSNITYRTINDYGALTPKMTGVME </pre>
<b>预测分子量</b>	25 kDa including tags
<b>氨基酸</b>	37 to 241
<b>标签</b>	His tag N-Terminus
<b>描述</b>	重组 <i>E. coli</i> fimC蛋白

#### 技术指标

Our **Abpromise guarantee** covers the use of **ab104677** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>应用</b>	SDS-PAGE Mass Spectrometry
<b>质谱法</b>	MALDI-TOF
<b>形式</b>	Liquid

## 制备和贮存

### 稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

## 常规信息


### 相关性

FimC is a member of the periplasmic chaperone family which functions in the chaperone-usher pathway and is indispensable in the biogenesis of the type 1 pilus fiber of *Escherichia coli*. It binds and interacts with fimH.

### 细胞定位

Periplasm

## 图片



15% SDS-PAGE analysis of 3µg ab104677.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors