

Recombinant *E. coli* glsA1 protein ab107143

1 图像

描述

| | | |
|-----------|--|--|
| 产品名称 | 重组 <i>E. coli</i> glsA1蛋白 | |
| 纯度 | > 95 % SDS-PAGE. Purified by using conventional chromatography techniques | |
| 表达系统 | <i>Escherichia coli</i> | |
| Accession | <u>P77454</u> | |
| 蛋白长度 | Full length protein | |
| 无动物成分 | No | |
| 性质 | Recombinant | |
| 种属 | <i>Escherichia coli</i> | |
| 序列 | | MLDANKLQQAVDQAYTQFHSLNGGQNADYIPFLANVPGQLAA VAIVTCDG NVYSAGDSDYRFALESISKVCTLALALEDVGPQAVQDKIGAD PTGLPFNS VIALELHGGKPLSPLVNAGAIATTSLINAENVEQRWQRILHI QQQLAGEQ VALSDEVNQSEQTTFHNRAIAWLLYSAGYLYCDAMEACDVY TRQCSTLL NTIELATLGATLAAGGVNPLTHKRVLQADNVPYILAEMMMEG LYGRSGDW AYRVGLPGKSGVGGGILAVVPGVMGIAAFSPPLDEDGNSVRG QKMVASVA KQLGYNVFKG |
| 预测分子量 | 35 kDa including tags | |
| 氨基酸 | 1 to 310 | |
| 标签 | His tag N-Terminus | |
| 描述 | 重组 <i>E. coli</i> glsA1蛋白 | |

技术指标

Our **Abpromise guarantee** covers the use of **ab107143** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----|-------------------------------|
| 应用 | Mass Spectrometry SDS-PAGE |
|----|-------------------------------|

| | |
|-----|-----------|
| 质谱法 | MALDI-TOF |
| 形式 | Liquid |

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00


Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息

序列相似性

Belongs to the glutaminase family.

图片



ab107143 at 3 µg analysed by 15% SDS PAGE.

SDS-PAGE - Recombinant *E. coli* glsA1 protein
(ab107143)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors