

Recombinant *E. coli* Multi Tag protein ab36791

1 References [1 图像](#)

描述

产品名称	重组 <i>E. coli</i> Multi Tag蛋白
生物活性	This protein lysate contains the commonly used epitope tags: HA, c-myc, Glu-Glu, GST, 6-His, HSV, T7, V5, VSV, E-tag, S-tag, FLAG (TM). These proteins were fused in the following order, GST, T7, HSV, c-myc, VSV, Glu-Glu, V5, E-tag, DYKDDDDK, S-tag and HA, 6xHis.
表达系统	<i>Escherichia coli</i>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	<i>Escherichia coli</i>
预测分子量	45 kDa
描述	重组 <i>E. coli</i> Multi Tag蛋白

技术指标

Our **Abpromise guarantee** covers the use of **ab36791** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
形式	Liquid
补充说明	<p>This protein lysate contains the commonly used epitope tags: HA, c-myc, Glu-Glu, GST, 6-His, HSV, T7, V5, VSV, E-tag, S-tag, FLAG (TM). These proteins were fused in the following order, GST, T7, HSV, c-myc, VSV, Glu-Glu, V5, E-tag, DYKDDDDK, S-tag and HA, 6xHis.</p> <p>This product is a bacterial lysate provided in 1 x reduced SDS-PAGE loading buffer. It is effective as a positive control for epitope tag antibodies in immunoblotting.</p> <p>For Western blotting, the suggested use of this reagent is 0.5 - 2.0 μl. It is recommended that the reagent be titrated for optimal performance for each application.</p>

制备和贮存

稳定性和存储	<p>Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.</p> <p>pH: 6.80</p>
--------	--

图片



All lanes : Recombinant *E. coli* Multi Tag protein (ab36791) at 0.5 μ l

Lane 1 : HA-HRP

Lane 2 : Myc-HRP

Lane 3 : Glu-Glu

Lane 4 : GST

Lane 5 : His

Lane 6 : HSV

Lane 7 : T7

Lane 8 : V5

Lane 9 : VSV

Lane 10 : S-Tag

Lane 11 : FLAG-HRP

Predicted band size: 45 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors