

# Recombinant Golden hamster GRP78 BiP protein ab92363

### 3 图像

#### 描述

---

产品名称	重组Golden仓鼠GRP78 BiP蛋白
纯度	> 80 % SDS-PAGE. Purified by Multi-step chromatography.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Hamster

#### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab92363** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot SDS-PAGE
形式	Liquid

#### 制备和贮存

---

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: Tris buffered saline, 0.77% (R*,R*)-1,4-Dimercaptobutan-2,3-diol, Glycerol (glycerin, glycerine), Sodium chloride
--------	--

#### 常规信息

---

功能	Probably plays a role in facilitating the assembly of multimeric protein complexes inside the endoplasmic reticulum. Involved in the correct folding of proteins and degradation of misfolded proteins via its interaction with DNAJC10, probably to facilitate the release of DNAJC10 from its substrate.
疾病相关	Autoantigen in rheumatoid arthritis.



**序列相似性**

Belongs to the heat shock protein 70 family.

**细胞定位**

Endoplasmic reticulum lumen. Melanosome. Cytoplasm. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

**图片**




**All lanes** : Anti-Recombinant GRP78 BiP polyclonal antibody

**Lane 1** : Molecular weight markers

**Lanes 2-3** : ab92363

Rabbit polyclonal to GRP78 BiP ([ab105019](#)) was tested with ab92363.



**All lanes :** Anti-KDEL (GRP78, GRP94) monoclonal antibody

**Lane 1 :** Molecular weight markers

**Lanes 2-3 :** ab92363

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors