

Recombinant Hepatitis B Virus E Antigen protein ab91273

6 References

描述

产品名称	重组乙型肝炎病毒E Antigen蛋白
生物活性	ab91273 is immunoreactive with sera of HBV-infected individuals. It contains the Hepatitis B Virus HBe adw immunodominant region.
纯度	> 95 % SDS-PAGE. Affinity purified.
表达系统	Escherichia coli
Accession	<u>Q17UT4</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
预测分子量	44 kDa including tags
氨基酸	8 to 165
标签	GST tag N-Terminus
额外的序列信息	The protein contains immunodominant region of HBV e antigen.

技术指标

Our **Abpromise guarantee** covers the use of **ab91273** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
	Western blot
	ELISA
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.79% Tris HCl, 50% Glycerol, 0.35% Sodium chloride, 0.307% Glutathione, 0.25% Sarkosyl
--------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

This product is an active protein and may elicit a biological response in vivo, handle with caution.

常规信息

相关性

Hepatitis B e-antigen (HBeAg) is a viral protein associated with HBV infections. Unlike the surface antigen, the e-antigen is found in the blood only when there are viruses also present. When the virus goes into "hiding," the e-antigen will no longer be present in the blood. HBeAg is often used as a marker of ability to spread the virus to other people (infectivity). Measurement of e-antigen may also be used to monitor the effectiveness of HBV treatment; successful treatment will usually eliminate HBeAg from the blood and lead to development of antibodies against e-antigen (anti-HBe). There are some types (strains) of HBV that do not make e-antigen; these are especially common in the Middle East and Asia. In areas where these strains of HBV are common, testing for HBeAg is not very useful.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors