

Recombinant Hepatitis E Virus ORF3 antigen protein ab49034

描述

产品名称	重组Hepatitis E Virus ORF3 antigen蛋白
纯度	> 95 % SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
序列	NPPDHSAPLGVTRPSAPPLPHVVDLPQLGPRR
氨基酸	92 to 123

技术指标

Our **Abpromise guarantee** covers the use of **ab49034** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

形式	Liquid
----	--------

制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle. pH: 7.2 Constituents: 48% Urea, 0.078% Beta mercaptoethanol, 0.316% Tris HCl
--------	---

常规信息

相关性	Hepatitis E virus (HEV) is a major human pathogen in the developing world. Very little is known about the basic biology of the virus. A small protein of unknown function, pORF3, is encoded by the third open reading frame of HEV.
-----	--

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors