

Recombinant HIV1 integrase protein ab49075

描述

产品名称	重组HIV1 integrase蛋白
生物活性	This protein reacts strongly with human HIV positive serum.
纯度	> 95 % SDS-PAGE.
表达系统	Escherichia coli
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant

技术指标

Our **Abpromise guarantee** covers the use of **ab49075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA Western blot
形式	Liquid
补充说明	This protein reacts strongly with human HIV positive serum. Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances. It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.2 Constituents: 0.2% Triton-X-100, 9% Urea, 0.395% Tris HCl, 50% Glycerol (glycerin, glycerine)
--------	--

常规信息

相关性	Integrase is an enzyme, a product of the HIV1 Gag/Pol gene; the other two enzymes being reverse
-----	---

transcriptase and protease. The human immunodeficiency virus (HIV-1) uses an enzyme, a so-called integrase, to carry out the integration of its viral DNA into the host chromosome thereby tricking the host cell machinery into making viral proteins. The HIV-1 integrase is a protein of 32 kDa and is composed of three domains. No cellular homologue of HIV integrase has been described, so potential inhibitors to this enzyme could be relatively nontoxic.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors