

# Recombinant human A-Raf protein ab79795

## 1 图像

### 描述

产品名称	重组人A-Raf蛋白
纯度	> 62 % SDS-PAGE. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
氨基酸	282 to 606

### 技术指标

Our **Abpromise guarantee** covers the use of **ab79795** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
形式	Liquid

### 制备和贮存

稳定性和存储	Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.0462% (R*,R*)-1,4-Dimercaptobutan-2,3-diol, 0.395% Tris HCl, 0.05% Tween, 50% Glycerol (glycerin, glycerine), 0.58% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--


### 常规信息

功能	Involved in the transduction of mitogenic signals from the cell membrane to the nucleus.
组织特异性	Predominantly in urogenital tissues.

## 序列相似性

Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. RAF subfamily.  
Contains 1 phorbol-ester/DAG-type zinc finger.  
Contains 1 protein kinase domain.  
Contains 1 RBD (Ras-binding) domain.

## 图片



10% SDS-PAGE showing ab79795 at approximately 63.6kDa (3μg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors