

Recombinant Human ABHD14B protein ab123190

1 图像

|           |  |  |
|-----------|--|--|
| 描述        |  |  |
| 产品名称      | 重组人ABHD14B蛋白   |  |
| 纯度        | > 95 % SDS-PAGE.<br>ab123190 is purified using conventional chromatography techniques.   |  |
| 表达系统      | Escherichia coli   |  |
| Accession | <u>Q96IU4</u>  |  |
| 蛋白长度      | Full length protein  |  |
| 无动物成分     | No   |  |
| 性质        | Recombinant  |  |
| 种属        | Human  |  |
| 序列        | MGSSHHHHHH SSGLVPRGSH MGSHEMAASVE<br>QREGTIQVQG QALFFREALP GSGQARFSVL<br>LLHGIRFSSE TWQNLGTLHR LAQAGYRAVA<br>IDLPGLGHSK EAAAPAPIGE LAPGSFLAAV<br>VDALELGPPV VISPSLSGMY SLPFLTAPGS<br>QLPGFVPVAP ICTDKINAAN YASVKTPALI<br>VYGDQDPMGQ TSFEHLKQLP NHRVLIMKGA<br>GHPCYLDKPE EWHTGLLDLFL QGLQ |  |
| 预测分子量     | 25 kDa including tags  |  |
| 氨基酸       | 1 to 210   |  |
| 标签        | His tag N-Terminus   |  |

技术指标

Our **Abpromise guarantee** covers the use of **ab123190** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用

Mass Spectrometry

SDS-PAGE

质谱法

MALDI-TOF

形式

Liquid

制备和贮存

## 稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 常规信息

### 功能

Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.

### 组织特异性

Ubiquitous. Detected in spleen, thymus, prostate, testis, ovary, small intestine, colon, peripheral blood leukocyte, heart, placenta, lung, liver, skeletal muscle, pancreas and kidney.


### 序列相似性

Belongs to the AB hydrolase superfamily. ABHD14 family.

### 细胞定位

Cytoplasm. Nucleus. Predominantly cytoplasmic.

## 图片



SDS-PAGE analysis of Recombinant Human ABHD14B protein (ab123190), under reducing conditions. Proteins visualized by coomassie blue stain.

Lane 1: Molecular Weight Standards

Lane 2: 3 µg Recombinant Human ABHD14B protein

SDS-PAGE - Recombinant Human ABHD14B protein (ab123190)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors