

Recombinant Human AG-3 protein ab108127

1 图像

描述	
产品名称	重组人AG-3蛋白
纯度	> 90 % SDS-PAGE. ab108127 was purified using conventional chromatography.
表达系统	Escherichia coli
Accession	<u>Q8TD06</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMGSMIAIKKEKRPPQTLSRG WGDDITWVQT YEEGLFYAQSKKPLMVIHHLEDCQYSQALKKVFAQNEEIQE MAQNKFIM LNLMHETTDKNLSPDGQYVPRIMFVDPSTVRADIAGRYSNR LYTYEPRD LPLLIENMKKALRLIQSEL
预测分子量	20 kDa including tags
氨基酸	22 to 166
标签	His tag N-Terminus
技术指标	
Our Abpromise guarantee covers the use of ab108127 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

常规信息


序列相似性

Belongs to the AGR family.

细胞定位

Secreted.

图片



SDS-PAGE analysis of Recombinant Human AG-3 protein (ab108127), under reducing conditions. Proteins visualized by coomassie blue stain.

Lane 1: Molecular Weight Standards

Lane 2: 3 µg Recombinant Human AG-3 protein

SDS-PAGE - Recombinant Human AG-3 protein
(ab108127)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors