

# Recombinant Human AHSA2 protein ab111452

## 描述

产品名称	重组人AHSA2蛋白
纯度	> 90 % SDS-PAGE.
表达系统	Escherichia coli
Accession	<b><u>Q719I0</u></b>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	38 kDa including tags
氨基酸	1 to 299
标签	His tag N-Terminus

## 技术指标

Our **Abpromise guarantee** covers the use of **ab111452** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
	SDS-PAGE
形式	Liquid

## 制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C.
	Constituents: 0.12% Tris buffer, 48.04% Urea, 1.1998% sodium dihydrogen phosphate

## 常规信息

功能	Co-chaperone that stimulates HSP90 ATPase activity.
序列相似性	Belongs to the AHA1 family.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors