

Recombinant Human ANAPC13 protein ab123152

1 图像


描述	
产品名称	重组人ANAPC13蛋白
纯度	> 90 % SDS-PAGE. ab123152 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u><b>Q9BS18</b></u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MASMTGGQQM GRGSHMDSEV QRDGRILDLI DDAWREDKLP YEDVAIPLNE LPEPEQDNGG TTESVKEQEM KWTDLALQYL HENVPPIGN
预测分子量	10 kDa including tags
氨基酸	1 to 74
标签	T7 tag N-Terminus
技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab123152</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 常规信息

功能	Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle. The APC/C complex acts by mediating ubiquitination and subsequent degradation of target proteins: it mainly mediates the formation of 'Lys-11'-linked polyubiquitin chains and, to a lower extent, the formation of 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains.
通路	Protein modification; protein ubiquitination.
序列相似性	Belongs to the APC13 family.
细胞定位	Nucleus.

## 图片



15% SDS-PAGE showing ab123152 at approximately 10kDa (3μg). (Molecular weight on SDS-PAGE will appear higher)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors