

Recombinant human ANGPTL7 protein ab109158

描述

产品名称	重组人ANGPTL7蛋白
生物活性	Activates ERK phosphorylation in THP1 cells.
纯度	> 90 % SDS-PAGE. ab109158 is an 0.2 µm filtered solution.
内毒素水平	< 0.100 Eu/µg
表达系统	HEK 293 cells
Accession	<u>O43827</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	45 kDa including tags
氨基酸	1 to 346
标签	DDDDK tag C-Terminus

技术指标

Our **Abpromise guarantee** covers the use of **ab109158** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	---

常规信息

组织特异性	Highly and specifically expressed in the cornea where is confined to the stromal layer.
序列相似性	Contains 1 fibrinogen C-terminal domain.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors