

Recombinant Human Annexin-13/ANXA13 protein ab105616

1 图像

描述	
产品名称	重组人Annexin-13/ANXA13蛋白
纯度	> 90 % SDS-PAGE. ab105616 was purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<b><u>P27216</u></b>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<b>MGSSHHHHHSSGLVPRGSHMGNRHAKASSPQGFVDVRDA KKLNKACKGM GTNEAAIIEILSGRTSDERQQIKQKYKATYGKEEVLKSEL SGNFEKTA LALLDRPSEYAARQLQKAMKGLGTDESVLIEVLCTR TNKEII AIKEAYQR LFDRSLES DVKGD TSGNLKKILV SLLQANRNEGDDVDKDLA G QDAKDLYD AGEGRWGTD ELAFNEVLAKRSYKQLRATFQAYQILIGKDIEE AIEEETSG DLQKAYLT LVRCAQDCEDYFAERLYKSMKGAGTDEETLIRIV VTRAEVDL QGIKAKFQEKYQKSLSDMVRSDTSGDFRKL LVAL LH</b>
预测分子量	38 kDa including tags
氨基酸	1 to 316
标签	His tag N-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab105616</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry

质谱法	MALDI-TOF
形式	Liquid
补充说明	Previously labelled as Annexin A13. This product was previously labelled as ANXA13


## 制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 0.0292% EDTA, 30% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
--------	---

## 常规信息

组织特异性	Gut specific.
序列相似性	Belongs to the annexin family. Contains 4 annexin repeats.
结构域	A pair of annexin repeats may form one binding site for calcium and phospholipid.
细胞定位	Cell membrane. Associated with the plasma membrane of undifferentiated, proliferating crypt epithelial cells as well as differentiated villus enterocytes.

## 图片



15% SDS-PAGE analysis of 3µg ab105616.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors