

Recombinant Human Apo-M protein ab93935

1 图像

描述	
产品名称	重组人Apo-M蛋白
纯度	> 90 % SDS-PAGE. ab93935 was purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>O95445</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MCPEHSQLTT LGVDGKEFPE VHLGQWYFIA GAAPTKEELA TFDPVDNIVF NMAAGSAPMQ LHLRATIRMK DGLCVPRKWI YHLTEGSTDL RTEGRPDMKT ELFSSSCP GG IMLNETGQGY QRFLLYNRSP HPPEKCV EEF KSLTSCLD SK AFLLT PRNQE ACELSNN
预测分子量	21 kDa
氨基酸	23 to 188
额外的序列信息	The letters in bold correspond to the His tag sequence.
技术指标	
Our <u>Abpromise guarantee</u> covers the use of ab93935 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

常规信息

功能	Probably involved in lipid transport. Can bind sphingosine-1-phosphate, myristic acid, palmitic acid and stearic acid, retinol, all-trans-retinoic acid and 9-cis-retinoic acid.
组织特异性	Plasma protein. Expressed in liver and kidney.
序列相似性	Belongs to the calycin superfamily. Lipocalin family. Highly divergent.
细胞定位	Secreted. Present in high density lipoprotein (HDL) and to a lesser extent in triglyceride-rich lipoproteins (TGRLP) and low density lipoproteins.

图片



15% SDS-PAGE analysis of 3µg ab93935

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors