

Recombinant Human ARL5A protein ab101943

1 图像

描述	
产品名称	重组人ARL5A蛋白
纯度	> 90 % SDS-PAGE. ab101943 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q9Y689</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMGSHMGILFTRIWRLFNHQE HKVIIVGLDN AGKTTILYQFSMNEVVHTSPTIGSNVEEIVINNTRFLMWDIG GQESLRSS WNTYYTNTEFVIVVVDSTDRERISVTREELYKMLAHEDLRKA GLLIFANK QDVKECMTVAEISQFLKLTSIKDHQWHIQACCALTGEGLCQG LEWMSRL KIR
预测分子量	23 kDa including tags
氨基酸	1 to 179
标签	His tag N-Terminus

技术指标	
Our Abpromise guarantee covers the use of ab101943 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.00348% PMSF, 0.0308% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine)

常规信息


相关性

ARL5A lacks ADP ribosylation enhancing activity. The protein encoded by this gene belongs to the ARF family of GTP binding proteins. With its distinctive nuclear/nucleolar localization and interaction with HP1alpha, the protein is developmentally regulated and may play a role(s) in nuclear dynamics and/or signaling cascades during embryonic development. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene has multiple pseudogenes.

细胞定位

Intracellular

图片



15% SDS-PAGE analysis of 3µg ab101943.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors