

Recombinant Human BAALC protein ab95394

1 图像

描述

产品名称	重组人BAALC蛋白
纯度	> 95 % SDS-PAGE. ab95394 is purified using conventional chromatography techniques.
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MCGGSRADA IEPRYYESWT RETESTWLTY TSDAPPSAA APDSGPEAGG LHSGMLEDGL PSNGVPRSTA PGGIPNPEKK TNCETQCPNP QSLSSGPLTQ KQNGLQTTEA KRDAKRMPAK EVTINVTDSI QQMDRSRRIT KNCVN
氨基酸	1 to 145

技术指标

Our **Abpromise guarantee** covers the use of **ab95394** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Mass Spectrometry SDS-PAGE
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 0.0584% EDTA, 20% Glycerol (glycerin, glycerine)
--------	---

常规信息

功能	May play a synaptic role at the postsynaptic lipid rafts by interacting with CAMK2A.
组织特异性	Expressed by neural and hematopoietic cells. Expression is found both in normal early progenitor cells and in the most immature type of blasts in acute leukemia but not in mature hematopoietic cells. Isoform 2 and isoform 6 are expressed in the brain. Isoform 2 shows a low expression in neuroectoderm-derived tissues such as adrenal gland and no expression in bone marrow, peripheral blood lymphocytes or lymph nodes, or in tumors or cancer cell lines of nonneural tissue origin.
翻译后修饰	Palmitoylation and myristoylation target the protein to the lipid rafts.
细胞定位	Cytoplasm. Membrane. Postsynaptic lipid rafts.

图片



15% SDS-PAGE showing ab95394 at approximately 17.7kDa (3 μ g).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors