

Recombinant Human BCMA protein ab139210

1 图像

描述	
产品名称	重组人BCMA蛋白
纯度	> 90 % SDS-PAGE. ab139210 was purified by using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q02223</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MGSRKINSEP LKDEFKNTGS GLLGMANIDL EKSRTGDEII LPRGLEYTVE ECTCEDCIKS KPKVDS DHCF PLPAMEEGAT ILVTTKTNDY CKSLPAALSA TEIEKSISAR
预测分子量	14 kDa including tags
氨基酸	78 to 184
标签	His tag N-Terminus
技术指标	
Our Abpromise guarantee covers the use of ab139210 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Mass Spectrometry SDS-PAGE
质谱法	MALDI-TOF
形式	Liquid
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol, 0.88% Sodium chloride

常规信息

功能	Receptor for TNFSF13B/BLYS/BAFF and TNFSF13/APRIL. Promotes B-cell survival and plays a role in the regulation of humoral immunity. Activates NF-kappa-B and JNK.
组织特异性	Expressed in mature B-cells, but not in T-cells or monocytes.
疾病相关	Note=A chromosomal aberration involving TNFRSF17 is found in a form of T-cell acute lymphoblastic leukemia (T-ALL). Translocation t(4;16)(q26;p13) with IL2.
序列相似性	Contains 1 TNFR-Cys repeat.
细胞定位	Cell membrane. Endomembrane system. Perinuclear Golgi-like structures.

图片



15% SDS-PAGE analysis of ab139210 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors