

Recombinant Human Blood Group Kell Antigen protein ab116891

1 图像

描述

产品名称	重组人Blood Group Kell Antigen蛋白
表达系统	Wheat germ
Accession	<u>P23276</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	<pre> MEGGDQSEEEPRERSQAGGMGLWSQESTPEERLPVEGSRPW AVARRVLT AILILGLLLCFSVLLFYNFQNCGPRPCETSVCLDLRDHYLAS GNTSVAPC TDFFSFACGRAKETNNSFQELATKNKNRLRRILEVQNSWHPG SGEEKAFQ FYNSCMDTLAIEAAGTGPLRQVIEELGGWRISGKWTSLNFR TLRLLMSQ YGHFPFFRAYLGPHPASPHTPVIQIDQPEFDVPLKQDQEKI YAQIFREY LTYLNQLGTLLGGDPSKVQEHSSLSISITSRLFQFLRPLEQR RAQGKLFQ MVTIDQLKEMAPAIDWLSCLQATFTPMSLSPSLSLVVHDVEY LKNMSQLV EEMLLKQRDFLQSHMILGLVVTLSPALDSQFQEARRKLSQKL RELTEQPP MPARPRWMKCVEETGTFEPTLAALFVREAFGPSTRSAAMKL FTAIRDAL ITRLRNLPWMNEETQNMAQDKVAQLQVEMGASEWALKPELAR QEYNDIQL GSSFQSVLSCVRSRLRARIVQSFLQPHPQHRWKVSPWDVNAY YSVSDHVV VFPAGLLQPPFFHPGYPRAVNFGAAGSIMAHELLHIFYQLLL PGGCLACD NHALQEHLCLKRHAAAFPLPSRTSFNDSLTFLENAADVGG AIALQAYS KRLLRHGETVLPDLSPQQIFFRSYAQVMCRKPSPQDSDH THSPHLR VHGPLSSTPAFARYFRCARGALLNPSSRCQLW </pre>

预测分子量 107 kDa including tags
氨基酸 1 to 732

技术指标

Our **Abpromise guarantee** covers the use of **ab116891** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 SDS-PAGE
Western blot
ELISA

形式 Liquid

制备和贮存

稳定性和存储 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.3% Glutathione, 0.79% Tris HCl


常规信息

功能 Zinc endopeptidase with endothelin-3-converting enzyme activity.

序列相似性 Belongs to the peptidase M13 family.

细胞定位 Cell membrane. Spans the erythrocyte membrane, and is attached to the underlying cytoskeleton.

图片



12.5% SDS-PAGE showing ab116891 at approximately 106.59kDa stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Blood Group Kell
Antigen protein (ab116891)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors