

Recombinant Human C8G protein ab113121

1 图像

描述		
产品名称	重组人C8G蛋白	
纯度	> 95 % SDS-PAGE. ab113121 was purified using conventional chromatography.	
表达系统	Escherichia coli	
Accession	<b><u>P07360</u></b>	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MGSSHHHHHHSSGLVPRGSHMQKPQRPRRPASPISTIQPKAN FDAQQFAG TWLLVAVGSACRFLQEQGHRAEATTLHVAPQGTAMAVSTFRK LDGICWQV RQLYGDTGVLGRFLLQARGARGAVHVVVAETDYQSFAVLYLE RAGQLSVK LYARSLPVSDSVLSGFQERVQEHLTEDQIFYFPKYGFCEAA DQFHLDE VRR	
预测分子量	23 kDa including tags	
氨基酸	21 to 202	
标签	His tag N-Terminus	

技术指标		
Our <b>Abpromise guarantee</b> covers the use of <b>ab113121</b> in the following tested applications.		
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.		
应用	Mass Spectrometry SDS-PAGE	
质谱法	MALDI-TOF	
形式	Liquid	

制备和贮存

<b>稳定性和存储</b>	<p>Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.</p> <p>pH: 8.00</p> <p>Constituents: 0.03% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.88% Sodium chloride</p>
<b>常规信息</b>	
<b>功能</b>	<p>C8 is a constituent of the membrane attack complex. C8 binds to the C5B-7 complex, forming the C5B-8 complex. C5-B8 binds C9 and acts as a catalyst in the polymerization of C9. The gamma subunit seems to be able to bind retinol.</p>
<b>序列相似性</b>	<p>Belongs to the calycin superfamily. Lipocalin family.</p>
<b>细胞定位</b>	<p>Secreted.</p>

## 图片



15% SDS-PAGE analysis of ab113121 (3μg)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors