

Recombinant human CCL1 protein ab9854

描述

产品名称	重组人CCL1蛋白
生物活性	Determined by its ability to chemoattract human monocytes using a concentration range of 10.0-30.0 ng/mL.
纯度	>= 98 % SDS-PAGE. = 98% by SDS-PAGE gel and HPLC analyses.
内毒素水平	< 1.000 Eu/μg
表达系统	Escherichia coli
Accession	<u>P22362</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	SKSMQVPFSR CCFSFAEQEI PLRAILCYRN TSSICSNEGL IFKLRGKEA CALDTVGWVQ RHRKMLRHCP SKRK
预测分子量	9 kDa

技术指标

Our **Abpromise guarantee** covers the use of **ab9854** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
形式	Lyophilized
补充说明	The biological activity of this product was determined by its ability to chemoattract total human T cell population using a concentration range of 10.0-100.0 ng/ml. This product was previously labelled as I-309

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
	This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	Reconstitute in water to a concentration of 0.1mg/ml.

常规信息

功能	Cytokine that is chemotactic for monocytes but not for neutrophils. Binds to CCR8.
序列相似性	Belongs to the intercrine beta (chemokine CC) family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors