

Recombinant human CCN3 protein ab50076

1 图像

描述	
产品名称	重组人CCN3蛋白
生物活性	Determined by a cell proliferation assay using BALB/c 3T3 cells. The expected ED ₅₀ for this effect is 1.0-2.0 µg/ml
纯度	> 95 % SDS-PAGE. ab50076 purity was determined also by HPLC. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MQVAATQRCP PQCPRGCPAT PPTCAPGVRA VLDGCSCCLV CARQRGESCS DLEPCDESSG LYCDRSADPS NQTGICTAVE GDNCVFDGVI YRSGEKFQPS CKFQCTCRDG QIGCVPRCQL DVLLPEPNCP APRKVEVPGE CCEKWICGPD EEDSLGGLTL AAYRPEATLG VEVSDSSVNC IEQTTEWTAC SKSCGMGFST RVTNRNRQCE MLKQTRLCMV RPCEQEPEQP TDKKGKKCLR TKKSLKAIHL QFKNCTSLHT YKPRFCGVCS DGRCTPHNT KTIQAEFQCS PGQIVKKPVM VIGTCTCHTN CPKNNEAFLQ EELKTTTRGK M
氨基酸	28 to 357

技术指标	
Our Abpromise guarantee covers the use of ab50076 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Western blot Functional Studies SDS-PAGE

形式 Lyophilized

制备和贮存

稳定性和存储

Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.

This product is an active protein and may elicit a biological response in vivo, handle with caution.

复溶

Reconstitute to 1mg/ml with distilled water.

常规信息

功能

Immediate-early protein likely to play a role in cell growth regulation.

组织特异性

Expressed in bone marrow, thymic cells and nephroblastoma. Increased expression in Wilms tumor of the stromal type.

序列相似性

Belongs to the CCN family.

Contains 1 CTCK (C-terminal cystine knot-like) domain.

Contains 1 IGFBP N-terminal domain.


Contains 1 TSP type-1 domain.

Contains 1 VWFC domain.

细胞定位

Secreted.

图片



Anti-CCN3 antibody (**ab10888**) at 1 µg/ml + Recombinant human CCN3 protein (ab50076) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (**ab97080**) at 1/5000 dilution

Exposure time: 1 minute

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors