

Recombinant human CD137 protein ab50088

描述	
产品名称	重组人CD137蛋白
生物活性	Determined by its inhibitory effect of the TNFSF9 (CD137 ligand) mediated stimulation of IL-8 production by human PBMC. About 90% of inibition was seen using a concentration of 1µg for both TNFSF9 and CD137.
纯度	> 98 % SDS-PAGE. The purity is greater than 98% by SDS-PAGE gel and HPLC analyses.
表达系统	Escherichia coli
Accession	<u><b>Q07011</b></u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MERTSLQDP CSNCPAGTFC DNNRNQICSP CPPNSFSSAG GQRTCDICRQ CKGVFTRKE CSSTSNAECD CTPGFHCLGA GCSMCEQDCK QGQELTKKGC KNCCFGTFND QKRGICRPWT NCSLDGKSVL VNGTKERDVV CGPSPADLSP GASSVTPPAP AREPGHS
预测分子量	18 kDa
氨基酸	19 to 184
技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab50088</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Inhibition Assay
形式	Lyophilized
制备和贮存	
稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -

20°C long term.

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 复溶

Reconstitute to 1mg/ml with distilled water.

## 常规信息

---

功能	Receptor for TNFSF14/4-1BBL. Possibly active during T cell activation.
组织特异性	Expressed on the surface of activated T-cells.
序列相似性	Contains 4 TNFR-Cys repeats.
细胞定位	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors