

# Recombinant Human CD30-L protein ab119813

## 描述

产品名称	重组人CD30-L蛋白
生物活性	Determined by its ability to stimulate Human IL8 production by Human PBMC using a concentration range of 10.0-25.0 ng/ml. Note: Results may vary with PBMC donors.
纯度	> 95 % SDS-PAGE. Greater than 95% by SDS-PAGE gel and HPLC analyses.
内毒素水平	< 0.100 Eu/μg
表达系统	CHO cells
Accession	<b><u>P32971</u></b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	HHHHHHHHPS PGGSGGQRTD SIPNSPDNVP LKGGNCSEDL LCILKRAPFK KSWAYLQVAK HLNKTKLSWN KDGILHGVRY QDGNLVIQFP GLYFIICQLQ FLVQCPNNSV DLKLELLINK HIKKQALVTV CESGMQTKHV YQNLSQFLLD YLQVNTTISV NVDTFQYIDT STFPLENVLS IFLYSNSD
氨基酸	63 to 234
标签	His tag N-Terminus

## 技术指标

Our **Abpromise guarantee** covers the use of **ab119813** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE HPLC
形式	Lyophilized
补充说明	Determined by its ability to stimulate Human IL8 production by Human PBMC using a concentration range of 10.0-25.0 ng/ml. Note: Results may vary with PBMC donors. Previously labelled as CD153.

## 制备和贮存

---

### 稳定性和存储

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Constituent: PBS

### 复溶

For lot specific reconstitution information please contact our Scientific Support Team.

## 常规信息

---

### 功能

Cytokine that binds to TNFRSF8/CD30. Induces proliferation of T-cells.

### 序列相似性

Belongs to the tumor necrosis factor family.

### 细胞定位

Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors