

Recombinant Human CD38 protein ab114253

1 图像

描述		
产品名称	重组人CD38蛋白	
表达系统	Wheat germ	
Accession	AAH07964	
蛋白长度	Full length protein	
无动物成分	No	
性质	Recombinant	
种属	Human	
序列	MANCEFSPVSGDKPCCRLSRAQLCLGVSILVLILVVVLAVV VPRWRQQW SGPGTTKRFPETVLARCVKYTEIHPEMRHVDCQSVWDAFKGA FISKHPCN ITEEDYQPLMKLGTQTVPCNKILLWSRIKDLAQFTQVQRDM FTLEDTLL GYLADDLTWCGEFNTSKINYQSCPDWRKDCSNNPVSFVKTV SRRFAEAA CDVVHAMLNGSRSKIFDKNSTFGSVEVHNLQPEKVQTLQAWV IHGGREDS RDLCDPTIKELSIISKRNIFSCCKNIYRDPKFLQCVKNPE DSSCTSEI	
预测分子量	59 kDa including tags	
氨基酸	1 to 300	

技术指标	
Our Abpromise guarantee covers the use of ab114253 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Western blot ELISA
形式	Liquid
补充说明	

制备和贮存

稳定性和存储

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

常规信息

功能

Synthesizes cyclic ADP-ribose, a second messenger for glucose-induced insulin secretion. Also has cADPr hydrolase activity. Also moonlights as a receptor in cells of the immune system.

组织特异性

Expressed at high levels in pancreas, liver, kidney, brain, testis, ovary, placenta, malignant lymphoma and neuroblastoma.

序列相似性

Belongs to the ADP-ribosyl cyclase family.


发展阶段

Preferentially expressed at both early and late stages of the B and T-cell maturation. It is also detected on erythroid and myeloid progenitors in bone marrow, where the level of surface expression was shown to decrease during differentiation of blast-forming unit E to colony-forming unit E.

细胞定位

Membrane.

图片



ab114253 analysed by 12.5% SDS-PAGE and stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors