

Recombinant Human CD99 protein (Fc Chimera His Tag) ab222999

1 图像

描述	
产品名称	重组人CD99蛋白(Fc Chimera His Tag)
纯度	> 95 % SDS-PAGE. Affinity purified
内毒素水平	< 1.000 Eu/μg
表达系统	Baculovirus infected insect cells
Accession	<u>P14209</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	ADPDGGFDLSDALPDNENKKPTAIPKKPSAGDDFDLGDVVD GENDDPRP PNPPKMPNPNPNHPSSSGSFSDADLADGVSGGEGKGGSDGG GSHRKEGE EADLEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMI SRTPEVTC VVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVV SVLTVLHQ DWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPP SRDELTKN QVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSGGS FFLYSKLT VDKSRWQQGNVFSCSVMHEALHNHYTQKSLSLSPGKHHHHHH
预测分子量	37 kDa including tags
氨基酸	23 to 122
标签	His tag C-Terminus
额外的序列信息	Extracellular domain fused to hlgG-His-tag at C-terminus.

技术指标

Our **Abpromise guarantee** covers the use of **ab222999** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 SDS-PAGE

形式 Liquid

制备和贮存

稳定性和存储 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.40

Constituents: 10% Glycerol (glycerin, glycerine), PBS

常规信息


功能 Involved in T-cell adhesion processes. It is involved in spontaneous rosette formation with erythrocytes.

序列相似性 Belongs to the CD99 family.

翻译后修饰 Extensively O-glycosylated.

细胞定位 Membrane.

图片



3 µg ab222999 was analysed by 15% SDS-PAGE.

SDS-PAGE - Recombinant Human CD99 protein (Fc Chimera His Tag) (ab222999)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors