

Recombinant Human CDKN2A/p16INK4a protein ab84075

[2 References](#) [1 图像](#)

描述

产品名称	重组人CDKN2A/p16INK4a蛋白		
纯度	> 95 % SDS-PAGE. Purity is greater than 95% by SDS-PAGE gel and HPLC analyses.		
内毒素水平	< 0.100 Eu/μg		
表达系统	Escherichia coli		
蛋白长度	Full length protein		
无动物成分	No		
性质	Recombinant		
种属	Human		
序列	MEPAAGSSME	PSADWLATAA	ARGRVEEVRA
	LLEAGALPNA	PNSYGRRPIQ	VMMMG SARVA
	ELLLLHGAEP	NCADPATLTR	PVHDAAREGF
	LDTLVVLHRA	GARLDVRDAW	GRLPVDLAE
	LGHRDVARYL	RAAAGGTRGS	NHARIDAAEG PSDIPD
预测分子量	17 kDa		
氨基酸	1 to 156		

技术指标

Our **Abpromise guarantee** covers the use of **ab84075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilized
补充说明	Endotoxin level is less than 0.1 ng per μg (1EU/μg).

制备和贮存


稳定性和存储	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. Constituents: 0.00268% PBS, 0.0154% DTT
复溶	Reconstitute with 1x PBS, pH 8.0. to a concentration of 0.1-1.0 mg/ml.

常规信息

细胞定位

Cytoplasmic and Nuclear

图片



Standard Curve for CDKN2A/p16INK4a (Analyte:

CDKN2A/p16INK4a protein (Human) (ab84075)); dilution range

1pg/ml to 1µg/ml using Capture Antibody **Mouse monoclonal**

[DCS50.1] to CDKN2A/p16INK4a (ab16123) at 1µg/ml and

Detector Antibody **Rabbit polyclonal to CDKN2A/p16INK4a**

(ab7962) at 0.5µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors