

Recombinant Human Chemerin protein ab104356

1 References

描述	
产品名称	重组人Chemerin蛋白
生物活性	Determined by its ability to chemoattract human immature dendritic cells using a concentration range of 1.0-100.0 ng/ml.
纯度	> 98 % SDS-PAGE. Purity > 98% by SDS-PAGE gel and HPLC analyses. Sterile filtered through a 0.2 micron filter.
内毒素水平	< 0.100 Eu/μg
表达系统	Escherichia coli
Accession	<u>Q99969</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	ELTEAQRRLGLQVALEEFHKHPPVQWAFQETSVESAVDTPFPA GIFVRLEF KLQQTSCRKRDWKKPECKVRPNRKRKCLACIKLGSEDKVLG RLVHCPIE TQVLREAEHQETQCLRVQRAGEDPHSFYFPGQFA
预测分子量	16 kDa
氨基酸	21 to 155

技术指标	
Our <u>Abpromise guarantee</u> covers the use of ab104356 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Functional Studies SDS-PAGE
形式	Lyophilized
补充说明	Determined by its ability to chemoattract human immature dendritic cells using a concentration range of 1.0-100.0 ng/ml. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. Do not vortex. This solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a

buffer containing a carrier protein (example 0.1% BSA) and store in working aliquots at -20°C to -80°C.

制备和贮存

稳定性和存储

Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Please see notes section.

复溶

Reconstitute in water to a concentration of 0.1-1.0 mg/ml. Do not vortex. This solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing a carrier protein (example 0.1% BSA) and store in working aliquots at -20°C to -80°C.

常规信息

组织特异性

Highly expressed in skin (basal and suprabasal layers of the epidermis, hair follicles and endothelial cells). Also found in pancreas, liver, spleen, prostate, ovary, small intestine and colon.

细胞定位

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors