

Recombinant human Chemerin protein ab124314

★★★★★ [1 Abreviews](#) [1 References](#)

描述

产品名称	重组人Chemerin蛋白
生物活性	ab124314 has full biological activity when compared to standards. The ED ₅₀ , as determined by its ability to chemoattract Human Chem23R transfected BaF3 mouse pro-B cells, is typically 4-20 ng/mL.
纯度	>= 98 % SDS-PAGE. Purity by SDS-PAGE: =98% Purity by HPLC: =98%
表达系统	Escherichia coli
Accession	Q99969
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MELTEAQRRG LQVALEEFHK HPPVQWAFQE TSVESAAVDTP FPAGIFVRLE FKLQQTSRK RDWKKPECKV RPNGRKRKCL ACIKLGSEDK VLGRLVHCPI ETQVLREAEE HQETQCLRVQ RAGEDPHSFY FPGQFAFS
预测分子量	16 kDa
氨基酸	21 to 157

技术指标

Our **Abpromise guarantee** covers the use of **ab124314** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE HPLC Functional Studies
形式	Lyophilized
补充说明	This product is manufactured by BioVision, an Abcam company and was previously called 4002 Chemerin, human recombinant. 4002-25 is the same size as the 25 µg size of ab124314.

Endotoxin level: <0.01 ng/μg.

制备和贮存

稳定性和存储

Shipped at 4°C. Store at -80°C. Avoid freeze / thaw cycle. For long term storage it is recommended to add a carrier protein on reconstitution (0.1% HSA or BSA).

Constituent: 0.2% Trifluoroacetic acid

This product is an active protein and may elicit a biological response in vivo, handle with caution.

复溶

Centrifuge vial before opening. When reconstituting the product, gently pipette and wash down the sides of the vial to ensure full recovery of the protein into solution. It is recommended to reconstitute the lyophilized product with sterile water to a concentration of 0.1mg/ml, which can be further diluted into other aqueous solutions.

常规信息

组织特异性

Highly expressed in skin (basal and suprabasal layers of the epidermis, hair follicles and endothelial cells). Also found in pancreas, liver, spleen, prostate, ovary, small intestine and colon.

细胞定位

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors