

Recombinant Human Chromogranin B protein ab93930

1 图像

描述

产品名称	重组人Chromogranin B蛋白
纯度	> 80 % SDS-PAGE. Purified using conventional chromatography techniques
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
氨基酸	21 to 677

技术指标

Our **Abpromise guarantee** covers the use of **ab93930** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
--------	---

常规信息


功能	Secretogranin-1 is a neuroendocrine secretory granule protein, which may be the precursor for other biologically active peptides.
组织特异性	Expressed in the adrenal medulla, and in pheochromocytoma. Not expressed in liver.
序列相似性	Belongs to the chromogranin/secretogranin protein family.

翻译后修饰

Extensively processed by limited proteolysis at conserved basic residues. Alternative processing are seen in different tissues.

细胞定位

Secreted. Neuroendocrine and endocrine secretory granules.

图片

15% SDS-PAGE showing ab93930 (3 µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors