

Recombinant Human CSRP2/Dlp1 protein ab126663

1 图像

描述	
产品名称	重组人CSRP2/Dlp1蛋白
纯度	> 85 % SDS-PAGE. ab126663 was purified using conventional chromatography techniques.
表达系统	Escherichia coli
Accession	<u>Q16527</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHH SSGLVPRGSH MGSHPVWGG GNKCGACGRT VYHAEVQCD GRSFHRCCFL CMVCRKNLDS TTVAIHDEEI YCKSCYGKKY GPKGYGYGQG AGTLNMDRGE RLGKIPESVQ PHRPTTNPNT SKFAQKYGGA EKCSRCGDSV YAAEKIIGAG KPWHKNCFRC AKCGKSLEST TLTEKEGEIY CKGCYAKNFG PKGFGYGQGA GALVHAQ
预测分子量	24 kDa including tags
氨基酸	1 to 193
标签	His tag N-Terminus

技术指标	
Our <b>Abpromise guarantee</b> covers the use of <b>ab126663</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid
补充说明	This product was previously labelled as CSRP2

## 制备和贮存

### 稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

## 常规信息

### 功能

Drastically down-regulated in response to PDGF-BB or cell injury, that promote smooth muscle cell proliferation and dedifferentiation. Seems to play a role in the development of the embryonic vascular system.

### 组织特异性

Highly expressed in the aorta, but not in heart and skeletal muscle.


### 序列相似性

Contains 2 LIM zinc-binding domains.

### 细胞定位

Nucleus.

## 图片



3ug by SDS-PAGE under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human CSRP2/Dlp1 protein (ab126663)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors