

Recombinant human CX3CL1 protein ab9820

描述

产品名称	重组人CX3CL1蛋白
生物活性	Determined by its ability to chemoattract human T cells using a concentration range of 5.0-10.0 ng/ml
纯度	> 97 % SDS-PAGE. Sterile filtered Greater than 97% pure by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g).
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	42 kDa

技术指标

Our **Abpromise guarantee** covers the use of **ab9820** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
形式	Lyophilized
补充说明	The biological activity of this product is determined by its ability to chemoattract human T-lymphocytes using a concentration range of 5.0-10.0 ng/ml.

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
	This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	For lot specific reconstitution information please contact our Scientific Support Team.

功能	The soluble form is chemotactic for T-cells and monocytes, but not for neutrophils. The membrane-bound form promotes adhesion of those leukocytes to endothelial cells. May play a role in regulating leukocyte adhesion and migration processes at the endothelium. Binds to CX3CR1.
组织特异性	Small intestine, colon, testis, prostate, heart, brain, lung, skeletal muscle, kidney and pancreas.
序列相似性	Belongs to the intercrine delta family.
翻译后修饰	A soluble short 95 kDa form may be released by proteolytic cleavage from the long membrane-anchored form. O-glycosylated with core 1 or possibly core 8 glycans.
细胞定位	Secreted and Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors