

### Recombinant human CXCL5 protein ab50039

#### 描述

产品名称	重组人CXCL5蛋白
生物活性	Biological Activity : Determined by its ability to chemoattract human neutrophils using a concentration range of 10.0-100.0 ng/ml.
纯度	> 98 % SDS-PAGE. Greater than 98% by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g).
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	LRELRCVCLQ TTQGVHPKMI SNLQVFAIGP QCSKVEVVAS LKNGKEICLD PEAPFLKKVI QKILDGGNKE N

#### 技术指标

Our **Abpromise guarantee** covers the use of **ab50039** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
形式	Lyophilized

#### 制备和贮存

稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.  This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	For lot specific reconstitution information please contact our Scientific Support Team.

#### 常规信息

功能	Involved in neutrophil activation. In vitro, ENA-78(8-78) and ENA-78(9-78) show a threefold higher chemotactic activity for neutrophil granulocytes.
序列相似性	Belongs to the intercrine alpha (chemokine CxC) family.
翻译后修饰	N-terminal processed forms ENA-78(8-78) and ENA-78(9-78) are produced by proteolytic cleavage after secretion from peripheral blood monocytes.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors