

Recombinant human Cyclin E1 + CDK3 protein (Active)
ab196131

2 图像

描述	
产品名称	重组人Cyclin E1 + CDK3蛋白(Active)
生物活性	Specific activity of ab196131 was determined to be 30 pmol/min/μg. Assay Conditions: Kinase buffer containing 1 mM DTT using Histone H1 (0.1 mg/mL) substrate and 20 μM ATP at 30°C for 50 min.
纯度	> 90 % SDS-PAGE. Affinity purified.
表达系统	Baculovirus infected Sf9 cells
Accession	<u>Q00526</u> <u>P24864</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
氨基酸序列 1	
种属	Human
预测分子量	60 kDa
标签	GST tag N-Terminus
额外的序列信息	Recombinant full length human Cdk3 (NM_001258)
氨基酸序列 2	
种属	Human
预测分子量	73 kDa
额外的序列信息	Recombinant full length Cyclin E1 (NM_001238)

技术指标	
Our Abpromise guarantee covers the use of ab196131 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE Functional Studies

形式	Liquid
补充说明	This is a complex of Cdk3 and Cyclin E1 Co-expressed in a Baculovirus-infected Sf9 cell expression system.


制备和贮存

稳定性和存储	<p>Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle.</p> <p>pH: 8.00</p> <p>Constituents: 0.63% Tris HCl, 0.64% Sodium chloride, 0.02% Potassium chloride, 0.12% Glutathione, 0.0462% (R*,R*)-1,4-Dimercaptobutan-2,3-diol</p> <p>This product is an active protein and may elicit a biological response in vivo, handle with caution.</p>
--------	---


常规信息

细胞定位	Cyclin E1: Nucleus.
------	---------------------

图片



Specific activity of ab196131 was determined to be 30 pmol/min/μg.



10% SDS-PAGE analysis of ab196131 stained with Coomassie Blue.

Lane 1: Cdk3 Lower band and Cyclin E1 upper band

Lane 2: Protein Marker

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors