

# Recombinant Human DcR1 protein ab151632

### 描述

---

产品名称	重组人DcR1蛋白
纯度	> 95 % SDS-PAGE. ab151632 is greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.
内毒素水平	< 1.000 Eu/μg
表达系统	HEK 293 cells
Accession	<b><u>O14798</u></b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	ATTARQEEVPQQTIVAPQQQRHSFKGEECPAGSHRSEHTGACN PCTEGVDY TNASNNEPSCFPCTVCKSDQKHKSSCTMTRDTCQCKEGTFR NENSPENC RKCSRCPSGEVQVSNCTSWDDIQCVEEFGANATVETPAAEET MNTSPGTP APAAEETMNTSPGTPAPAAEETMTTSPGTPAPAAEETMTTSP GTPAVDHH HHHH
预测分子量	22 kDa including tags
氨基酸	26 to 221
标签	His tag C-Terminus

### 技术指标

---

Our **Abpromise guarantee** covers the use of **ab151632** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE HPLC
形式	Lyophilized

### 制备和贮存

---

## 稳定性和存储

Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.

pH: 7.2

Constituents: 64% Sodium chloride, 3% Sodium phosphate monobasic, monohydrate, 28% disodium;hydrogen phosphate;dodecahydrate

## 复溶

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## 常规信息

---

### 功能

Receptor for the cytotoxic ligand TRAIL. Lacks a cytoplasmic death domain and hence is not capable of inducing apoptosis. May protect cells against TRAIL mediated apoptosis by competing with TRAIL-R1 and R2 for binding to the ligand.

### 组织特异性

Higher expression in normal tissues than in tumor cell lines. Highly expressed in peripheral blood lymphocytes, spleen, skeletal muscle, placenta, lung and heart.

### 序列相似性

Contains 3 TNFR-Cys repeats.

### 翻译后修饰

N-glycosylated and O-glycosylated.

### 细胞定位

Cell membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors