

Recombinant Human DcR2 protein (His tag) ab220541

1 图像

描述	
产品名称	重组人DcR2蛋白(His tag)
纯度	> 92 % SDS-PAGE.
内毒素水平	< 1.000 Eu/μg
表达系统	HEK 293 cells
Accession	<b><u>Q9UBN6</u></b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	ATIPRQDEVPPQTVAPQQRRSLKEEECPAGSHRSEYTGACN PCTEGVDY TIASNNLPSCLLCTVCKSGQTNKSSCTTTTRDTVCQCEKGSFQ DKNSPEMC RTCRTGCPRGMVKVSNCTPRSDIKCKNESAASSTGKTPAAEE TVTTILGM LASPYH
预测分子量	18 kDa
氨基酸	56 to 211
标签	His tag C-Terminus
额外的序列信息	Accession # NP_003831.

技术指标	
Our <b><u>Abpromise guarantee</u></b> covers the use of <b>ab220541</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	SDS-PAGE
形式	Lyophilized

制备和贮存	
稳定性和存储	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 7.40  
Constituents: 95% PBS, 5% Trehalose

Lyophilized from 0.22 µm filtered solution.

#### 复溶

Reconstitute with sterile deionized water to a concentration of 200 µg/ml.

#### 常规信息

##### 功能

Receptor for the cytotoxic ligand TRAIL. Contains a truncated death domain and hence is not capable of inducing apoptosis but protects against TRAIL-mediated apoptosis. Reports are contradictory with regards to its ability to induce the NF-kappa-B pathway. According to PubMed:9382840, it cannot but according to PubMed:9430226, it can induce the NF-kappa-B pathway.

##### 组织特异性

Widely expressed, in particular in fetal kidney, lung and liver, and in adult testis and liver. Also expressed in peripheral blood leukocytes, colon and small intestine, ovary, prostate, thymus, spleen, pancreas, kidney, lung, placenta and heart.


##### 序列相似性

Contains 1 death domain.  
Contains 3 TNFR-Cys repeats.

##### 细胞定位

Membrane.

#### 图片



SDS-PAGE analysis of ab220541 under reducing conditions, stained overnight with Coomassie Blue. As a result of glycosylation, the reduced protein migrates as 33-40 kDa.

SDS-PAGE - Recombinant Human DcR2 protein  
(His tag) (ab220541)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors