

Recombinant Human DcR3 protein ab157279

1 图像

描述	
产品名称	重组人DcR3蛋白
生物活性	ab157279 binds Human LIGHT. Does not bind Human FasL.
纯度	> 95 % SDS-PAGE.
内毒素水平	< 0.100 Eu/μg
表达系统	HEK 293 cells
Accession	<b><u>O95407</u></b>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Human
序列	MRALEGPGLS LLCLVLALPA LLPVPAVRGV AETPTYPWD AETGERLVCA QCPPGTQVQR PCRRDSPTTC GPCPPRHYTQ FWNYLERCY CNVLCGEREE EARACHATHN RACRCRTGFF AHAGFCLEHA SCPPGAGVIA PGTPSQNTQC QPCPPGTQSA SSSSSEQCQP HRNCTALGLA LNVPGSSSHD TLCTSCTGFP LSTRVPGAEE CERAVIDFVA FQDISIKRLQ RLLQALEAPE GWGPTP
预测分子量	58 kDa including tags
氨基酸	1 to 246

技术指标	
Our <b><u>Abpromise guarantee</u></b> covers the use of <b>ab157279</b> in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
应用	Functional Studies SDS-PAGE
形式	Lyophilized
补充说明	ab157279 binds Human LIGHT. Does not bind Human FasL.

## 制备和贮存

### 稳定性和存储

Shipped at 4°C. Store at -20°C.

Constituent: 99% PBS

### 复溶

Reconstitute with 50µl sterile water to a concentration of 1 mg/ml. After reconstitution, prepare aliquots and store at -20°C. Avoid freeze/thaw cycles.

## 常规信息

### 功能

Decoy receptor for the cytotoxic ligands TNFSF14/LIGHT and TNFSF6/FASL. Protects against apoptosis.

### 组织特异性

Detected in fetal lung, brain and liver. Detected in adult stomach, spinal cord, lymph node, trachea, spleen, colon and lung. Highly expressed in several primary tumors from colon, stomach, rectum, esophagus and in SW480 colon carcinoma cells.


### 序列相似性

Contains 4 TNFR-Cys repeats.

### 细胞定位

Secreted.

## 图片



ELISA analysis of ab157279.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors