

Product datasheet

Recombinant Human Density Regulated Protein ab111640

1 图像

描述

产品名称	重组人Density Regulated蛋白
纯度	> 85 % SDS-PAGE. ab111640 is purified using conventional chromatography.
表达系统	Escherichia coli
Accession	<u>O43583</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MGSSHHHHHSSGLVPRGSHMAADISESSGADCKGDPRNS AKLDADYPLR VLYCGVCSLPTEYCEYMPDVAKCRQWLEKNFPNEFAKLTVEN SPKQEAGI SEGQGTAGEEEEEKKQKRGGRGQIKQKKKTVPQKVTIKIPR AKKKYVTR VCGLATFEIDLKEAQRFFAQKFSCGASVTGEDEIIIQGDFTD DIIDVIQE KWPEVDDDSIEDLGEVKK
预测分子量	24 kDa including tags
氨基酸	1 to 198
标签	His tag N-Terminus
描述	重组人Density Regulated蛋白

技术指标

Our **Abpromise guarantee** covers the use of **ab111640** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
形式	Liquid

制备和贮存

稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.04% DTT, 20% Glycerol (glycerin, glycerine), 0.88% Sodium chloride

Contains 89% Tris-HCl buffer.

常规信息

功能

May be involved in the translation of target mRNAs by scanning and recognition of the initiation codon. Plays a role in the modulation of the translational profile of a subset of cancer-related mRNAs when recruited to the translational initiation complex by the oncogene MCTS1.

组织特异性


Highly expressed in heart and skeletal muscle and moderately expressed in the brain, placenta, liver and pancreas. Weakly expressed in the lung and kidney.

序列相似性

Belongs to the DENR family.

Contains 1 SUI1 domain.

图片



SDS-PAGE of ab111640 (3µg) under reducing condition and visualized by coomassie blue stain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors