

Recombinant Human DFNA5/GSDME protein ab116824

1 图像

描述	
产品名称	重组人DFNA5/GSDME蛋白
表达系统	Wheat germ
Accession	<u>O60443</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
序列	MFAKATRNFLREVDADGDLIAVSNLNDSDKLQLLSLVTKKKR FWCWQRPK YQFLSLTLGDVLIEDQFSPVWVESDFVKYEGKFANHVSGTL ETALGKVK LNLGGSSRVESQSSFGLTRKQEVDLQQLIRDSAERTINLRNP VLQQVLEG RNEVLCVLTKITTMQKCVISEHMQVEEKCGGIVGIQTKTVQ VSATEDGN VTKDSNVVLEIPAATTIAYGVIELYVKLDGQFEFCLLRGKQG GFENKKRI DSVYLDPLVFREFAFIDMPDAAHGISSQDGPLSVLKQATLLL ERNFHPFA ELPEPQQTALSDIFQAVLFDDELLMVLEPVCDDLVSGLSPTV AVLGELKP RQQQDLVAFLQLVGCSLQGGCPGPEDAGSKQLFMTAYFLVSA LAEMPDSA AALLGTCCKLQIIPTLCHLLRALSDDGVSLEDPTLTPLKDT ERFGIVQR LFASADISLERLKSSVKAVILKDSKVFPLLLCITLNGLCALG REHS
预测分子量	81 kDa including tags
氨基酸	1 to 496
标签	GST tag N-Terminus
技术指标	

Our **Abpromise guarantee** covers the use of **ab116824** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	SDS-PAGE
	Western blot
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
	pH: 8.00
	Constituents: 0.3% Glutathione, 0.79% Tris HCl

常规信息

组织特异性	Expressed in cochlea. Low level of expression in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas, with highest expression in placenta.
疾病相关	Defects in DFNA5 are the cause of deafness autosomal dominant type 5 (DFNA5) [MIM:600994]. DFNA5 is a form of sensorineural hearing loss. Sensorineural deafness results from damage to the neural receptors of the inner ear, the nerve pathways to the brain, or the area of the brain that receives sound information.
序列相似性	Belongs to the gasdermin family.

图片



12.5% SDS-PAGE stained with Coomassie Blue showing ab116824 at approximately 80.63 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors