

Recombinant Human DNAJA1 protein ab113174

2 图像

描述

产品名称	重组人DNAJA1蛋白
纯度	> 90 % SDS-PAGE. ab113174 was purified by multi-step chromatography.
表达系统	Escherichia coli
Accession	<u>P31689</u>
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
预测分子量	45 kDa
氨基酸	1 to 394

技术指标

Our **Abpromise guarantee** covers the use of **ab113174** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot SDS-PAGE
形式	Liquid
补充说明	This product was previously labelled as HDJ2


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: 99% Tris buffered saline, 0.02% (R*,R*)-1,4-Dimercaptobutan-2,3-diol
--------	---

常规信息

功能	Co-chaperone of Hsc70. Seems to play a role in protein import into mitochondria.
序列相似性	Contains 1 CR-type zinc finger. Contains 1 J domain.
细胞定位	Membrane.

图片



SDS-PAGE - Recombinant Human DNAJA1 protein (ab113174)

SDS-PAGE analysis


Lane 1: MW marker

Lane 2: 0.5 µg DNAJA1 protein (ab113174)

Lane 3: 1 µg DNAJA1 protein (ab113174)

Lane 4: 2 µg DNAJA1 protein (ab113174)

Lane 5: 5 µg DNAJA1 protein (ab113174).



Western blot - Recombinant Human DNAJA1 protein (ab113174)

Western Blot analysis

Lane 1: MW marker

Lane 2: DNAJA1 (ab113174) probed with anti-DNAJA1 monoclonal antibody.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors